


Installation Guide for the TJ Transfer Case Lowering Kit

Part #001900100 (TJ 97-02, Skin Packed) Part #001900190 (TJ 97-02, Boxed)

> Part #001900103 (TJ 03-06, Skin Packed) Part #001900193 (TJ 03-06, Boxed)

- 1. 1997-2002 TJ Wranglers- Install the six spacers (1.25"OD x 1.25" tall) between the frame and the cross member with the flat end of the spacer up against the frame. The 1/2" lock washer goes onto the 1/2" x 3" bolt first followed by the tapered washer with the flat end of the tapered washer against the lockwasher.
- 2. Put the tires and wheels back on if needed and lower the TJ to the ground.



- 1. 2003-current TJ Wranglers- Install the six spacers (1.25"0D x 1.25" tall) between the frame and the cross member with the flat end of the spacer up against the frame. The 1/2" lock washer goes onto the M12 x 1.75 x 70 bolt first followed by the 1/2" flat washer.
- 2. Put the tires and wheels back on if needed and lower the TJ to the ground.